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
Fig. 1 Surface area of Gram-negative bacteria

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
Fig. 2 Principle of endotoxin measurement using Lysate reagent (encircled in blue)

(a) Coagulation (gel-clot method): If a clot is apparent on inversion of the tube after 37 °C for one hour, the result is positive. (b) Coagulation (turbidimetric method): The amount of gel resulted from the cleavage of coagulogen is optically measured. (c) Coloration (colorimetric method): pNA produced from the cleavage of substrate by clotting enzyme is optically measured. (d) Luminescence (bioluminescence method): D-luciferin produced from the cleavage of substrate by clotting enzyme is measured.

Table 1 Comparison of reaction time and detection limit in endotoxin detection methods⁵⁾

Method	Reaction Time /min	Detection Limit /EU mL ⁻¹
End-point turbidimetric	15	0.01
End-point chromogenic	15	0.1
End-point chromogenic (Long reaction time)	30	0.01
Turbidimetric kinetic	138	0.0005
Kinetic chromogenic	90	0.002
Bioluminescence (This method)	15	0.0005

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
Fig. 3 The three-dimensional structure of North American firefly luciferase and amino acid substitution positions of high-luminescence mutants³⁾


Fig. 4 Time-course changes of bioluminescence intensities of mutant luciferase FM and wild-type luciferase⁴⁾

Reactions contained 100 μL of luciferase (0.25 mg mL^{-1}) and 100 μL of substrate (10^{-8} mol ATP and 10^{-9} mol D-luciferine), both in tris/HCl buffer.

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
Fig. 5 Comparison of endotoxin detection curves by bioluminescence and chromogenic methods⁵⁾

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
Fig. 6 Effects of various concentrations of sodium chloride, potassium chloride, calcium chloride, ammonium chloride, and sodium acetate on the activity of the wild-type luciferase⁷⁾

The relative luminescence intensity (%) was calculated by dividing the recorded intensity of luminescence in the presence of an indicated salt by that in the absence of a salt. The values are represented as means $\pm$ SD ($n = 3$).


Fig. 7 Effect of NaCl concentration on luminescence intensities of luciferase FM

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
Fig. 8 Effect of NaCl on bioluminescence method using luciferase FM⁸⁾

Luminescence intensities were measured at various concentration of luciferin.


Fig. 9 Luminescence intensities of bioluminescence method using luciferase FM under dialysis condition⁸⁾

Luminescence intensities were measured at various concentration of endotoxin.

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Table 2 The effect of Val288 substitutions on luciferase activity in the presence of 140 mmol L⁻¹ sodium chloride⁷⁾

Luciferase	Substituted amino acid	Molecular weight /g mol ⁻¹	Polarity and charge	Relative luminescence intensity/ %
Wild type	Val	117.15	Nonpolar	44.0 ± 0.6
Mutant	Gly	75.07	Nonpolar	16.3 ± 0.1
Mutant	Ala	89.09	Nonpolar	21.9 ± 1.3
Mutant	Pro	115.13	Nonpolar	19.6 ± 0.4
Mutant	Ile	131.17	Nonpolar	67.2 ± 2.9
Mutant	Leu	131.17	Nonpolar	46.8 ± 0.6
Mutant	Met	149.21	Nonpolar	47.3 ± 0.7
Mutant	Phe	165.19	Nonpolar	46.2 ± 0.8
Mutant	Trp	204.23	Nonpolar	28.9 ± 0.9
Mutant	Ser	105.09	Polar	16.3 ± 0.1
Mutant	Thr	119.12	Polar	25.6 ± 0.5
Mutant	Cys	121.16	Polar	46.5 ± 1.3
Mutant	Asn	132.12	Polar	29.3 ± 0.5
Mutant	Gln	146.15	Polar	32.8 ± 0.3
Mutant	Tyr	181.19	Polar	38.3 ± 0.9
Mutant	Lys	146.19	Positive	23.8 ± 1.8
Mutant	His	155.15	Positive	31.1 ± 1.4
Mutant	Arg	174.20	Positive	40.8 ± 0.8
Mutant	Asp	133.10	Negative	16.1 ± 0.4
Mutant	Glu	147.13	Negative	25.1 ± 1.0

Table 3 The effect of Glu488 substitutions on luciferase activity in the presence of 140 mmol L⁻¹ sodium chloride⁷⁾

Luciferase	Substituted amino acid	Molecular weight /g mol ⁻¹	Polarity and charge	Relative luminescence intensity/ %
Wild type	Val	147.13	Negative	44.0 ± 0.6
Mutant	Gly	75.07	Nonpolar	50.2 ± 0.2
Mutant	Ala	89.09	Nonpolar	43.4 ± 1.1
Mutant	Pro	115.13	Nonpolar	54.5 ± 1.2
Mutant	Ile	117.15	Nonpolar	78.7 ± 1.3
Mutant	Leu	131.17	Nonpolar	58.0 ± 0.5
Mutant	Met	131.17	Nonpolar	62.0 ± 0.5
Mutant	Phe	149.21	Nonpolar	73.1 ± 0.7
Mutant	Trp	165.19	Nonpolar	56.6 ± 1.1
Mutant	Ser	204.23	Nonpolar	40.0 ± 1.2
Mutant	Thr	105.09	Polar	50.9 ± 1.0
Mutant	Cys	119.12	Polar	66.8 ± 1.7
Mutant	Asn	121.16	Polar	69.5 ± 1.1
Mutant	Gln	132.12	Polar	30.9 ± 1.1
Mutant	Tyr	146.15	Polar	34.8 ± 1.4
Mutant	Lys	181.19	Polar	40.2 ± 0.7
Mutant	His	146.19	Positive	22.9 ± 1.3
Mutant	Arg	155.15	Positive	32.8 ± 0.5
Mutant	Asp	174.20	Positive	43.1 ± 1.0
Mutant	Glu	133.10	Negative	68.1 ± 1.2

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
Fig. 10 Inhibition of the activity of wild-type and mutant luciferases by sodium chloride⁷⁾

The luminescence intensity of luciferase reactions in the absence or presence of 140 mmol L^{-1} sodium chloride was measured. The values are represented as means $\pm$ SD ($n = 3$).

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
Fig. 12 The effect of sodium chloride concentration on the activity of wild-type luciferase, luciferase FM and luciferase CR⁷⁾


Fig. 11 The three-dimensional structure of *Photinus pyralis* luciferase and the locations of mutations in the luciferase CR⁷⁾

The active sites with a substrate analog (dehydroluciferylsulfamoyl adenosine) are indicated in the right box. The amino acid substitutions (red) in the luciferase CR are indicated. This figure was created based on the crystal structure (Protein Data Bank ID code 4G36) of *Photinus pyralis* luciferase.


Fig. 13 Effect of NaCl on bioluminescence method using luciferase⁸⁾

Luminescence intensities were measured at various concentrations of luciferin.


Fig. 14 Luminescence intensities of bioluminescence method using luciferase CR under dialysis condition⁸⁾

Luminescence intensities were measured at various concentrations of endotoxin.

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Fig. 15 Two-step reaction type bioluminescence ET analyzer


Fig. 16 Evaluation of the bioluminescence method by other groups⁹⁾

Six third parties measured endotoxin concentrations under dialysis and water conditions.

Fig. 17 One-step reaction type bioluminescence ET analyzer

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
Fig. 18 ET-detection biochemiluminescence methods: Two-step reaction versus One-step reaction


Fig. 19 Relationship of measured endotoxin in distilled water and endotoxin in dialysis solution by using bioluminescence endotoxin detection based on one-step addition of reagent

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
Fig. 20 Developments of bioluminescence ET detection method and analyzers

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Development of a Bioluminescence Detection Method for Endotoxin and Its Application to Dialysate Analysis

Satoshi YAWATA¹, Kenichi NODA¹, Ai SHIMOMURA¹, Atsumu ODA¹,
Satoshi ARAKAWA¹, Hiromitsu HACHIYA^{*1} and Akio KURODA²

^{*} E-mail : hachiya@toadkk.co.jp

¹ DKK-TOA CORPORATION, 613, Kita-Iriso, Sayama-shi, Saitama, 350-1388

² Graduate School of Integrated Sciences for Life, Hiroshima University, 1-3-1, Kagamiyama, Higashi-Hiroshima-shi, Hiroshima, 739-8530

(Received December 23, 2022; Accepted March 27, 2023)

A novel endotoxin detection method was developed that combines a bioluminescence method using a high-luminant mutant luciferase and a lysate reaction. This method reached the detection limit of 0.0005 EU mL⁻¹ in a 15-min duration. A salt-resistant mutant luciferase was developed whose sensitivity was not decreased by sodium chloride (NaCl) in dialysate; it was applied to the endotoxin analyzer for dialysate measurement. It was useful for dialysate measurement with a detection limit of 0.0003 EU mL⁻¹ at 20 min and exhibited other properties under dialysis conditions. A thermostable luciferase was also developed and applied to a single-reaction type endotoxin analyzer. The detection limit was 0.001 EU mL⁻¹ for 10 min, enabling a more rapid endotoxin detection. This novel endotoxin detection method has rapidly been becoming popular in the field of dialysate measurement because of its excellent features.

Keywords: bioluminescence; endotoxin; detection; luciferase; lysate; dialysate.